

AVANOS

DIGESTIVE HEALTH CORFLO*
Product Catalogue

**YOUR PREFERRED
ENTERAL PARTNER**



Key to Product Code

PRODUCT CATEGORY	ENFit® PRODUCT CODE
CORFLO® NG/NI Feeding Tubes	42-xxxx
CORTRAK® 2 EAS Feeding Tubes	42-xxxxTRAK2
CORFLO® For Endoscopic Placement	50-xxxx
FARRELL® Valve System	44 -4100
PEG Tubes and Accessories	50-xxxx
Jejunostomy (J-Tube) Extension Tubes (for use with CORFLO® PEG Tubes)	50-xxxx



**ENFit® Connectors
are available for all
feeding tube categories.**

**YOUR CLEAR CONNECTION
TO ENTERAL SAFETY**

CORTRAK* 2

Enteral Access System (EAS)

CORTRAK* 2 Enteral Access System (EAS)

PRODUCT CODE	DESCRIPTION	PACKAGED
20-0950	CORTRAK* 2 Enteral Access System (EAS) All-In-One Monitor, Integrated Visual Display Terminal (VDT), Touch Screen and Embedded Computing System	1 per case
	Smart Receiver Unit	Included
	Power Supply	Included
	Power Cord	Included
	Interconnect Cable	Included
	Operator's Manual with Service Guide	Included
	Tip Guide Booklet	Included
20-0951	CORTRAK* 2 EAS Accessory Kit	1 per case

CORTRAK* 2 Enteral Access System (EAS) Replacement Parts

PRODUCT CODE	DESCRIPTION	PACKAGED
20-0952	Smart Receiver Unit	1 per case
20-0966	Interconnect Cable	1 per case
20-0953	Li-ION Battery Pack	1 per case

CORTRAK* 2 Enteral Access System (EAS) Accessory Items

PRODUCT CODE	DESCRIPTION	PACKAGED
20-0911	Printer Paper Rolls (125 per Roll)	3 per case
20-0916	Printer Replacement Battery (QL220 Plus)	1 per case
20-0920	Stand (Includes SRU/RU Cover Holder)	1 per case
20-0928	QL220 Plus Printer (includes 2 rechargeable batteries)	1 per case
20-0929	QL220 Plus Quad Charger	1 per case
20-0930	Plastic Covers (Sleeves) for SRU/RU (4 rolls, 25 bags per roll)	100 per case
20-0955	Stand Bracket Kit	1 per case

CORTRAK* 2 Enteral Access System (EAS) Feeding Tubes

CORTRAK* 2 Nasogastric/Nasointestinal Feeding Tube with Electromagnetic Stylet

PRODUCT CODE	DESCRIPTION	PACKAGED
42-9368TRAK2	8 FR 36"/91 cm Length	10 tubes per case
42-9361TRAK2	10 FR 36"/91 cm Length	10 tubes per case
42-9362TRAK2	12 FR 36"/91 cm Length	10 Tubes per case
42-9438TRAK2	8 FR 43"/109 cm Length	10 tubes per case
42-9431TRAK2	10 FR 43"/109 cm Length	10 tubes per case
42-9558TRAK2	8 FR 55"/140 cm Length	10 tubes per case
42-9551TRAK2	10 FR 55"/140 cm Length	10 tubes per case



CORFLO*

NG/NI Feeding Tubes

All CORFLO* NG/NI Feeding Tubes are fully radiopaque, polyurethane dual port for the administration of medication and flushing, maintaining closed system, pre-inserted stylet, water activated lubricant on tube tip and internal lumen, regular cm markings, unique anti-clog tip

Ask for the **YELLOW** Tube...

CORFLO* Nasogastric/Nasointestinal Feeding Tube with Stylet

Anti-clog, bowed tip bolus, non-weighted, adult-pediatric-neonatal, with stylet

PRODUCT CODE	DESCRIPTION	PACKAGED
42-9156	6 FR 15"/38 cm Length	10 tubes per case
42-9225*	5 FR 22"/56 cm Length	10 tubes per case
42-9226	6 FR 22"/56 cm Length	10 tubes per case
42-9228	8 FR 22"/56 cm Length	10 tubes per case
42-9366	6 FR 36"/91 cm Length	10 tubes per case
42-9368	8 FR 36"/91 cm Length	10 tubes per case
42-9361	10 FR 36"/91 cm Length	10 tubes per case
42-9438	8 FR 43"/109 cm Length	10 tubes per case
42-9431	10 FR 43"/109 cm Length	10 tubes per case
42-9432	12 FR 43"/109 cm Length	10 tubes per case

CORFLO* Nasogastric/Nasointestinal Feeding Tube

Anti-clog, bowed tip bolus, non-weighted, adult-pediatric-neonatal, without stylet

PRODUCT CODE	DESCRIPTION	PACKAGED
42-1362	12 FR 36"/91 cm Length	10 tubes per case

CORFLO* Nasogastric/Nasointestinal Feeding Tube with Stylet

Anti-clog, smooth shape bolus, 3 gm* adult weighted, with stylet

PRODUCT CODE	DESCRIPTION	PACKAGED
42-7366	6 FR 36"/91 cm Length, 1.5 gm*	10 tubes per case
42-7368	8 FR 36"/91 cm Length	10 tubes per case
42-7361	10 FR 36"/91 cm Length	10 tubes per case
42-7558	8 FR 55"/140 cm Length	10 tubes per case

* Approximate weight

CORFLO* Nasogastric/Nasointestinal Feeding Tube with Stylet

Anti-clog, smooth shape bolus, pediatric weighted (1gm*), with stylet

PRODUCT CODE	DESCRIPTION	PACKAGED
42-8226	6 FR 22" /56 cm Length	10 tubes per case

CORFLO* Nasogastric/Nasointestinal Feeding Tube with Stylet

Pill shape bolus, 3 gm* adult weighted, with stylet

PRODUCT CODE	DESCRIPTION	PACKAGED
42-2438	8 FR 43" /109 cm Length	10 tubes per case
42-2431	10 FR 43" /109 cm Length	10 tubes per case
42-2551	10 FR 55" /140 cm Length	10 tubes per case

* Approximate weight

CORFLO* Nasointestinal Endoscopically Placed Feeding Tube

60" CORFLO* Anti-Clog Feeding Tube, 300 cm Teflon Coated Flexible Guide Wire (0.035" diameter) in dispenser sheath, and 10 FR, 22" length Nasal/Oral Transfer Tube

PRODUCT CODE	DESCRIPTION	PACKAGED
51-4602	CORFLO* ENDO/F NI 12 FR 7 gm* weighted	1 per case
51-9601	CORFLO* ENDO/F NI 10 FR Non-weighted	1 per case

* Approximate weight

CORFLO* Specialty Small Bowel Feeding (SBF) Tube Kits

Anti-clog, bowed tip bolus, non-weighted, with stylet

PRODUCT CODE	DESCRIPTION	PACKAGED
42-9431	10 FR 43" /109 cm Length	10 kits per case

CORFLO* Nasogastric/Nasointestinal Feeding Tube with Stylet

Anti-clog, bowed tip bolus, non-weighted adult tube, with stylet

PRODUCT CODE	DESCRIPTION	PACKAGED
42-9558	8 FR 55" /140 cm Length	10 tubes per case
42-9551	10 FR 55" /140 cm Length	10 tubes per case



FARRELL* Valve

Enteral Decompression System

FARRELL* Valve Closed Enteral Decompression System

Gastric pressure relief device, designed to help those who suffer from poor gastric motility, pain, and bloating. Maximises enteral feeding by providing a channel to constantly decompress the stomach, allowing the stomach to feed at its own pace.

PRODUCT CODE	DESCRIPTION	PACKAGED
44-4100	FARRELL* Valve FR 15"/38 cm Length Closed Enteral Decompression System	30 per case



Johnson Test Papers

CE marked for human gastric aspirate for the assessment of correct enteral tube placement

- Distinctive colours at key decision values - pH 4-6
- Tethered flip top lid minimises spillage and contamination
- Instructions for use directly on container
- Expiry date and lot number visible on container

PRODUCT CODE	DESCRIPTION	PACKAGED
103.3CE	pH Indicator Strip 0-6	100 strips per container
113.3CE	pH Indicator Strip 2-9	100 strips per container



CORGRIP*

NG/NI Tube Retention System

CORGRIP* System

for the securement of Nasogastric and Nasojejunal tubes

Kit contents:

- White Insertion Catheter with umbilical tape and magnet
- Yellow Retrieval Catheter with magnet
- SSL™ Clamp (Slot-Slide-Lock) with imprinted French size (2 clamps per kit)
- Two opening devices
- Lubricant

CORGRIP* NG/NI Tube Retention System

PRODUCT CODE	DESCRIPTION	PACKAGED
25-008	For use with 8 FR NG/NI feeding tubes	10 per case
25-010	For use with 10 FR NG/NI feeding tubes	10 per case
25-012	For use with 12 FR NG/NI feeding tubes	10 per case
25-014	For use with 14 FR NG/NI feeding tubes	10 per case
25-016	For use with 16 FR NG/NI feeding tubes	10 per case
25-018	For use with 18 FR NG/NI feeding tubes	10 per case



CORFLO*

PEG Kits

The CORFLO* PEG Kit is a kit for endoscopic placement of a gastrostomy tube, available in Pull and Push version (OTW Over the wire/Seldinger technique).

The kit contains the following components:

- Polyurethane transparent PEG tubing with a yellow radiopaque stripe and traction removable internal bolster
- Ring-shaped or conical internal retention bolster (see table)
- Two types of external retention disc
- Fast release clamp (Roberts clamp)
- Feeding connector
- Towel and drape
- Scalpel
- Jelly
- Stainless steel puncture needle and trocar Gauge 18
- Kit is ETO sterilised

The "E1-type" of kits have:

- An additional G14 puncture needle with teflon trocar
- No jelly

The "E2-type" of kits have:

- A blue coated pull wire
- 1 G18 puncture needle, except Fr 12 kit has a Gauge 14 puncture needle with teflon coated trocar
- No scalpel

The product item# 50-5020E2 contains in addition to E2-type of kits:

- 3 needles: Gauge 14, 18, 25
- Scalpel, syringe and gauzes

Pull Technique

PRODUCT CODE	DESCRIPTION	FR SIZE/BOLSTER DESIGN	PACKAGED
50-6020E1	CORFLO* PEG Kit	20 FR-Ring	2 per case
50-6016E1	CORFLO* PEG Kit	16 FR-Ring	2 per case
50-6012E1	CORFLO* PEG Kit	12 FR-Ring	2 per case
50-6020E2	CORFLO* PEG Kit	20 FR-Ring	2 per case
50-6016E2	CORFLO* PEG Kit	16 FR-Ring	2 per case
50-6012E2	CORFLO* PEG Kit	12 FR-Ring	2 per case
50-6520	CORFLO* PEG Kit	20 FR-Ring	2 per case
50-6516	CORFLO* PEG Kit	16 FR-Ring	2 per case
50-6512	CORFLO* PEG Kit	12 FR-Ring	2 per case
50-5520	CORFLO* PEG Kit	20 FR-Conical	2 per case
50-5020E2	CORFLO* PEG Kit	20 FR-Conical	2 per case

Push Technique

PRODUCT CODE	DESCRIPTION	FR SIZE/BOLSTER DESIGN	PACKAGED
50-4020E1	CORFLO* PEG Kit	20 FR-Ring	2 per case
50-4016E1	CORFLO* PEG Kit	16 FR-Ring	2 per case



CORFLO[®] PEG Feeding Connector Repair Kits

For Use With CORFLO[®] PEG Tubes

PRODUCT CODE	DESCRIPTION	PACKAGED
50-6120	CORFLO [®] PEG Y-Adapter Repair Kit–20 FR	5 per case
50-6116	CORFLO [®] PEG Y-Adapter Repair Kit–16 FR	5 per case
50-6112	CORFLO [®] PEG Y-Adapter Repair Kit–12 FR	5 per case

CORFLO[®] PEJ-Tubes

For Use With CORFLO[®] PEG Tubes

The CORFLO[®] PEJ is a jejunal extension to be used in combination with the CORFLO[®] PEG tube. Weighted Jejunal Tube with CORFLO[®] Anti-Clog Feeding Port. The tubing has centimeter markings for adjustment and placement confirmation and has a Grip-Tip to ease endoscopic placement, Adjust-a-Sleeve, Y-Adapter, Silicone Plug, and Stylet with Cap, Gastric access for aspiration, decompression and medication is maintained while using the CORFLO[®] PEJ.

PRODUCT CODE	DESCRIPTION	PACKAGED
50-7366	6 FR, 36"/91 cm; 3 gm ⁺ Weighted For Use With 16 FR PEG Tube	5 per case
50-7368	8 FR, 36"/91 cm; 3 gm ⁺ Weighted For Use With 20 FR PEG Tube	5 per case
50-7361	10 FR, 36"/91 cm; 3 gm ⁺ Weighted For Use With 20 FR PEG Tube	5 per case

⁺ Approximate weight

CORFLO[®] Over-the-Wire PEJ

For Use With 20 FR CORFLO[®] PEG Tubes

The CORFLO[®] Over-the-Wire PEJ is a jejunal extension to be used in combination with the CORFLO[®] PEG Fr 20 tube, allowing the Selindinger technique. Non-Weighted PEJ Tube with CORFLO[®] Anti-Clog Feeding Port and Adjust-a-Sleeve, Adapter, 260 cm Teflon Coated Flexible Guide Wire (0.035" diameter) in Dispenser Sheath and Silicone Plug.


PRODUCT CODE	DESCRIPTION	PACKAGED
50-1438	8 FR, 36"/91 cm; 3 gm ⁺ Weighted For Use With 20 FR PEG Tube	2 per case
50-1431	10 FR, 36"/91 cm; 3 gm ⁺ Weighted For Use With 20 FR PEG Tube	2 per case

AVANOS[®] Digestive Health - Your Preferred Enteral Partner

Unique enteral solutions including the FARRELL[®] gastric pressure relief system and the CORTRAK[®] Enteral Access System encourage clinicians to choose Avanos as their Preferred Enteral Partner.

Our research and development teams continuously work to bring innovative technologies to market, partnering with clinicians to create clinical and economic solutions that improve patient safety and outcomes.

Our products are supported by our respected value added services, including face to face training and e-learning through our Account Managers and Digestive Health education portal.



**From the UK, please call Avanos
at 0800 917 65 85 - fax 0800 169 02 35
or send an email to
customerservice.uk.ie@avanos.com**

**From N-Ireland - please call Sangers Surgical
at 028 9056 7135 - fax 028 9070 5623
or send an email to Hospital@sangersaah.co.uk**

**From IRL, please call Uniphar Hospital Services
at +353 1 468 84 56 - fax +353 1 466 07 82
or send an email to unipharhospitalservices@uniphar.ie**

**From all other countries, please call Avanos
at +32 2 700 68 51 - fax +32 2 711 26 91
or send an email to customerservice.export@avanos.com**

AVANOS

* Registered Trademark or Trademark of Avanos Medical, Inc., or its affiliates. © 2019 AVNS. All rights reserved. HC114-01-UK EtQ COPY-02165
ENFit® is a Registered Trademark or Trademark of Global Enteral Device Supplier Association, Inc. Used with permission.
The COLOR YELLOW is a Registered Trademark or Trademark of Avanos Medical, Inc. or its affiliates.
For Healthcare professional use only.